

## HOW TO USE INJECTION PUMP



<p><b>Connect the pump to electric power</b></p>
<p style="text-align: center;"><b>Load the syringe into the syringe slot</b></p> <p>Use the button at the back of the pump to move the pusher block to the left until it is pressing on the syringe barrel.</p> <p style="text-align: center;">See table below for explanation on how to move the pusher block</p>
<p style="text-align: center;"><b>Set up the syringe diameter:</b> Press SET DIAMETER</p> <p style="text-align: center;">See table below for most common syringe diameters</p>
<p style="text-align: center;"><b>Set REFILL rate to a maximum value</b></p> <p style="text-align: center;">This is the speed with which the pusher block moves backwards</p>
<p style="text-align: center;"><b>Set the TARGET volume in ml!</b></p> <p style="text-align: center;">This is the total volume of the contrast that will be injected.</p> <p>For a rat of 200 g, this will be 200 <math>\mu</math>l (dose 0.1 mmol/kg) or 400 <math>\mu</math>m (dose 0.2 mmol/kg)</p> <p style="text-align: center;"><b>NOTE: If this is the first injection, add 100 <math>\mu</math>l to the total volume to flush the dead space inside the cannula.</b></p>
<p style="text-align: center;"><b>Set the INFUSE rate in <math>\mu</math>m/min!</b></p> <p style="text-align: center;">To compute the infuse rate, use the formula below.</p> <p style="text-align: center;">In general, multiply target V (in <math>\mu</math>l): by 6 for a 10 sec or by 4 for a 15 sec injection</p>
<p style="text-align: center;"><b>Push RUN to start the infusion</b></p> <p style="text-align: center;">Make sure the display says INFUSE when you are ready to start</p> <p style="text-align: center;">If it says INTERRUPT, you will get a wrong volume</p>

REFILL = pusher block moves back  
 INFUSE = pusher block moves forward

Set Target Volume, Infuse Rate and Refill Rate by:  
 Press SET + TARGET VOLUME/INFUSE RATE/ REFILL RATE + DESIRED NUMBER + ENTER

Infuse rate ( $\mu$ l/min) =

$$\frac{\text{Target Volume (ml)} * 1000 \mu\text{l} / 1\text{ml}}{10 \text{ s} * 1 \text{ min} / 60 \text{ s}} = \frac{\text{Target Volume } (\mu\text{l})}{(1/6)}$$

1 ml syringe = 4.78 mm  
 3 ml Omnifix Braun = 9.68 mm  
 5 ml Omnifix = 12.45 mm  
 Insulin syringe (8mm tip) = 3.5 mm